

CellerateRX[®]Surgical



MEDCILLARY[®]

What is CellerateRX[®] Surgical?

- **CellerateRX[®] Surgical Activated Collagen[®] Powder** is a medical hydrolysate of Type I bovine collagen and contains no additives.
- **US FDA 510K cleared unclassified medical device**
- **Indicated** for the management of:
 - Surgical wounds
 - Traumatic wounds
 - Partial and full-thickness wounds
 - First and second degree burns
- **Sterile product and packaging** for use in a **sterile environment.**




CellerateRX[®] Surgical Powder

- No adverse events reported
- **Contraindications:** Known sensitivity to bovine collagen
- **Store at room temperature and protect from freezing**
- **3-year shelf life**



Supports the Body's Wound Healing Process

CellerateRX[®] Surgical aids in the management of surgical wounds



- Native Collagen is broken down into hydrolyzed collagen fragments
- CellerateRX[®] Surgical is approximately 1/100th the size of native collagen
- Hydrolyzed collagen fragments do not have to be broken down by the body before use, so the benefits of collagen are provided to the wound.

CellerateRX[®] Surgical Powder Illustration


Illustration of CellerateRX[®] Surgical Powder as small fragments



Illustration of CellerateRX[®] Surgical Powder + Moisture



Type I Collagen Hydrolyzed Fragments




CellerateRX[®] Surgical Is Appropriate For Many Surgical Specialties and Surgical Wounds

Orthopedic
Spine
Plastic/Reconstructive
Sports medicine
Trauma
General
Podiatric
Vascular
Obstetric/Gynecologic
Neurologic
Cardiovascular
Burn




Partial knee replacement



Open distal femur fracture



Open tibia and fibula fracture



Total abdominal hysterectomy



Hip replacement with comorbidities



Non-healing clavicle

CellerateRX® Surgical: Available in two delivery forms

Surgical Powder: Available in 1 gram and 5 gram vials

- Provides the optimum environment for healing; helps protect the wound from physical and bacterial injury (applied above fascia layer)
- Rapidly forms a gel in the presence of wound exudate
 - CellerateRX Surgical Powder 1 gram (List Price: \$895)
 - CellerateRX Surgical Powder 5 gram (List Price: \$2,895)



Surgical Gel: 28 gram tube

- Ideal for covering incisions (sutures/staples) at closing and during wound healing
 - CellerateRX Gel 28 gram tube (List Price: \$295)



How Is CellerateRX[®] Surgical Powder Applied?

1. Prepare the surgical wound site.
2. Apply CellerateRX[®] Surgical Powder to the surgical wound above the fascia layer, covering the entire wound bed.
3. Cover the surgical wound as per physician's protocol.

